CustomBone Service

Italian **IT**

English **EN**

French **FR**

Spanish **ES**

German **DE**

Portuguese **PT**

Russian **RU**

Slovak **SK**

Arabic **AR**

	<p>Non ristilizzare Do not Re-sterilize Ne pas Restériliser No volver a esterilizar Nicht erneut sterilisieren Não Re-esterilizar Не подвергать повторной стерилизации Nesterilizovať opätovne مەبەقەت قەدەخ زەنگەهەي ال</p>
	<p>Non riutilizzare Do not Re-use Ne pas Réutiliser No reutilizar Nicht wiederverwenden Não Reutilizar Не использовать повторно Nepoužiť opätovne مەدەخەتسەا قەدەخ زەنگەهەي ال</p>
	<p>Non utilizzare se la confezione è danneggiata Do Not Use if Package is Damaged Ne pas Utiliser si l'Enveloppe est Endommagée No utilizar si el envase está dañado Das Produkt nicht verwenden, wenn Beschädigungen der Verpackung festgestellt wurden Não Utilizar se a Embalagem estiver Danificada Не использовать при поврежденной упаковке Nepoužiť, ak je balenie poškodené قەسەلات تۆبەخەل تەنەك ادا دەخەتسەي ال</p>
	<p>Sterile, incluso il metodo di sterilizzazione (radiazioni ionizzanti) Indicates a medical device that has been sterilized using irradiation Indique un dispositif médical stérilisé par irradiation Esterilizado mediante radiaciones ionizantes Steril, einschließlich Strahlensterilisation Indica que o dispositivo médico foi esterilizado com radiações Стерильно, прошло стерилизацию, в т. ч. излучением Sterilné, vrátane sterilizačnej metódy (ionizovaná radiácia) (فەزەيۇمەتا تەخەخەلەل) مەبەقەتسەلا قەيخەرط نەبەزەتەم مەبەقەم</p>
	<p>Sterile, incluso il metodo di sterilizzazione (vapore o calore secco) Indicates a medical device that has been sterilized using steam or dry heat Indique un dispositif médical stérilisé par vapeur ou chaleur sèche Esterilizado con vapor o calor seco Steril, samt Sterilisationsmethode (Dampf oder Heißluft) Indica que o dispositivo médico foi esterilizado com vapor ou calor seco Стерильно, прошло стерилизацию паром или сухим жаром Sterilné, vrátane sterilizačnej metódy (para alebo suché teplo) (فەخەخەلەل تەخەخەلەل) مەبەقەتسەلا قەيخەرط نەبەزەتەم مەبەقەم</p>

	<p>Utilizzare entro Use-by Date Date d'expiration Fecha de vencimiento Verfallsdatum Data de Validade Использовать до Použiť do Použití do لإلاخ مادختسرالل</p>
	<p>Fabbricante Manufacturer Fabricant Fabricante Hersteller Fabricante Производитель Výrobca عزصنملا</p>
	<p>Numero di lotto Lot Number Numéro de Lot Número de lote Chargennummer Número do Lote Код партии Číslo výrobnéj dávky جياتسرالل عن حش</p>
	<p>Numero paziente Patient Number Numéro de Patient Número del paciente Patientennummer Número do Paciente Код пациента Číslo pacienta ضنملا جقدر</p>

	<p>Numero di catalogo Catalogue Number Numéro de Catalogue Número de catálogo Katalognummer Número de Catálogo Номер по каталогу Katalógové číslo چول اتشلا مؤر</p>
	<p>Numero di pezzi Number of Pieces Nombre de Pièces Número de piezas Stückzahl Número de Unidades Количество изделий Počet kusov عطقلا ددع</p>
	<p>Attenzione, vedere le istruzioni per l'uso Caution - See Instructions for Use Attention Atención: consultar las instrucciones de uso Siehe Gebrauchsanweisung Cuidado Внимание, См. инструкции по использованию Pozor: pozri NÁVOD NA POUŽITIE مادتختسالا تاداشرلا ةءارق بجي ،ءابتن</p>
	<p>Consultare le istruzioni operative Consult Instructions for Use Instructions d'Emploi Consultar las instrucciones de uso Bedienungsanleitung beachten Consultar as Instruções de Uso См. инструкции по применению Nahliadnuť do operačných inštrukcií لي غتختسالا تاداشرلا ةءارق بجي</p>

CUSTOMBONE SERVICE

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CUSTOMBONE SERVICE

CUSTOM-MADE IMPLANT FOR THE RECONSTRUCTION OF CRANIAL DEFECTS

INDICATIONS FOR USE

DESCRIPTION

Hydroxyapatite is contained in human bones in percentages ranging from 65% to 80% and it is one of the most important elements of human bone structure. CustomBone Service, patient specific implants for the reconstruction of cranial/craniofacial defects, are made of porous bio-mimetic hydroxyapatite with a chemical composition and structure that resembles the mineral component of human bones.

This biomaterial is highly porous, with trabecular structure and it is composed of pores with the following characteristics:

- ✓ macro-pores, which favor the migration and attachment of the cells responsible for bone regeneration;
- ✓ interconnecting pores, which are essential for ensuring that the cells colonising the implant receive sufficient nutrients and physiological liquids;
- ✓ micro-pores, which facilitate the absorption of physiological liquids.

This material is completely biocompatible and promotes rapid cell colonisation resulting in the formation of new bone.

The implants are designed and produced by Fin-Ceramica Faenza according to the surgeon's specifications and based on the patient's CT scan data, obtained through a standardized protocol.

During the pre-operative planning phase, the surgeon must approve the final design of the device.

For each device a specific conformity declaration is made available.

INTENDED USE

CustomBone Service is intended to replace bony voids in the cranial and/or craniofacial skeleton.

This device may be used both for adult and paediatric patients (for children above 2 years of age).

CustomBone Service implants are suitable for reconstruction of cranial defects resulting from:

- ✓ trauma and vascular pathologies, either associated or non-associated to cranial decompression;
- ✓ removal of tumours;
- ✓ reabsorption of autologous bone;
- ✓ rejection of other prosthetic materials;
- ✓ congenital malformations.

AVAILABLE FORMATS

CustomBone Service implants are patient specific: the devices are designed based upon the prescription of a qualified surgeon and related patient CTscan data. Each patient-specific implant is supplied with a back-up device. The implants are supplied in sterile packages.

INSTRUCTIONS FOR USE

These instructions are intended as guidelines for the use of CustomBone Service implants; they are not designed to replace or change the standard procedures for the treatment of cranial defects.

Clinical results of implants used for the reconstruction of cranial defects depend on several factors.

When choosing an implant and the surgical technique to be used, the following factors need to be considered: the patient's age and general clinical conditions, the bone quality, the possibility to achieve tight contact between the implant and the vascularized bone tissue, certainty of the complete filling of the defect as well as the possibility to obtain a correct and sufficient primary stabilisation of the implant.

Prior to the surgical intervention, make sure that the identification data on the implant match exactly with the patient ones as detailed in the patient's documents and medical records.

In addition, before the surgery, check all documents related to the project to make sure that the implant match exactly with the bone gap to be treated.

Pre-operative treatment

As per standard surgical practice, prior to the surgery the patient should be subjected to common antibiotic treatment. In patients that are allergic to specific antibiotics, an alternative treatment should be considered.

It is necessary to carefully verify that there is no infection or inflammation at the time of the operation.

Intra-operative aspects

Before applying the CustomBone Service implant, it is recommended to soak the product in an antibiotic solution in order to avoid any possible infection during the surgery.

Soaking in an antibiotic is a precautionary measure to prevent aerial contamination during the prosthesis affixing procedures while awaiting final fixation. Use of the antibiotic does not affect the final result. In cases of antibiotic allergies it is recommended to handle the prosthesis with sterile gauze.

Once the cranial defect has been exposed, it is necessary to remove any fibrotic or dural tissues from the bone edges to ensure maximum surface contact between the bone and the implant. The integration of the implant is highly favored when in contact with the greatest amount of vascularised bone tissue. Avoid exerting excessive pressure on the implant while positioning it; incorrect handling could in fact lead to the prosthesis damage. In order to stabilise the implant, prepare the suture holes on the edges of the bone specular to the ones on the edges of the prosthesis. Then, secure the implant using suture thread with a diameter less than 2 mm (non-reabsorbable). The central part of the implant also has holes for dural suspension; the decision to carry out this procedure is left to the discretion of the surgeon, and it depends on the overall clinical picture of the patient.

Once completed the implant fixation, the surgical site should be closed according to standard surgical procedures.

Post-operative procedure

In accordance with standard post-surgical procedures, an appropriate antibiotic therapy should be administered.

The surgeon must provide the patient with all the indications for a correct post-operative recovery, in relation to the localisation and entity of the defect as well as to the overall clinical picture of the patient.

The patient should be advised to pay particular attention to avoid direct traumas in the area of the implant during the first months after the surgery. If there are no post-operative complications and the post-operative indications given by the surgeon are correctly followed, the primary stabilization of the implant generally occurs after 6-8 weeks, through the continuing process of osseointegration that must be verified through a diagnostic evaluation.

CONTRAINDICATIONS

Use of CustomBone Service prostheses should be avoided in presence of inflammatory conditions and manifest infections of the surgical area to be treated.

ted, as well as in insufficiently vascularized sites. The CustomBone Service prostheses should not be implanted in case of dura mater lacerations. In addition, the device should not be used in patients suffering from a proven hypersensitivity to calcium phosphates or serious forms of demineralization diseases, in cases of chronic hypertension, hemocoagulation disorders, in the event of infections resulting from previous surgery or in patients treated with steroid or anti-coagulation therapies.

No side effects have been noted in the use of CustomBone Service where tumours have been present. Nonetheless, patients who undergo treatment with this device following the removal of a neoplasm should be monitored.

The product must not be implanted in areas in which skin tissue is not sufficient to cover the implant entirely.

WARNINGS

In pre-operative and intra-operative phases the device must be handled with the greatest care, avoiding any manoeuvres that might damage the device. If, during the surgery, the device is accidentally damaged or contaminated, complete the surgery using the back-up implant.

During the post-operative phase (6-8 weeks) any stress on the implant area should be avoided.

The device is designed and made to measure, exclusively for the patient indicated on the requesting physician's prescription.

Thus, it is recommended not to modify the custom-made device in any way. Any change made to the device supplied shall be the sole and exclusive responsibility of the surgeon. Any modification, however, should be made with extreme caution, with irrigation and the use of low speed diamond drills.

Indications and warnings given to the patient by the surgeon for the post-operative period are extremely important; in particular, patients should be warned to avoid direct traumas in the implant area. Severe impact on the implant area might lead to complications such as the mobilisation and/or fracture of the prosthesis.

In order to prevent any possible risk of bacterial infection, implants should be positioned in areas where a completely and totally covering with skin tissue is possible.

In order to promote a correct device stabilization, it is recommended to carefully evaluate - both during project validation and during surgery - any hypertension conditions or brain herniations which may hinder proper implant positioning.

Notably, in cases where reconstruction of the frontal bone might involve direct

exposure with the frontal sinuses open onto the nasal tracts, surgery is intended to keep the implant/airways completely separate.

As the device cannot withstand direct and indirect loads, once it is implanted it is not possible to apply skin expanders or use other similar techniques.

Do not use resorbable sutures for fixation.

The use of CustomBone Service implant is reserved exclusively to qualified medical staff.

CustomBone Service are single-use products; any residues must be disposed of. The product may not be re-sterilised.

The product should not be used if the packaging has been damaged. Do not use the product if it has expired.

Prior to surgery, it is imperative to check the labels on the implant and to make sure all the patient's labels are present and matching the traceability information on the implant and on all the packaging parts.

The package contains additional labels with the product's traceability details.

In the event of another surgery it is possible to use the device back-up only after evaluating the congruence between the gap and the device.

ADVERSE EFFECTS

No adverse effects have been reported following use of the product according to correct indications.

STERILIZATION

All CustomBone Service products are supplied sterile (i.e sterilized with gamma rays or steam). The products are single-use and may not be re-sterilised.

Its total or partial reuse may involve the risk of cross contamination and the danger of infection in the implant site.

STORAGE

The product must be stored in a cool and dry area, and should be protected from direct light and heat sources (+10 °C / +40 °C).

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FIN-CERAMICA FAENZA S.p.A.