




# Performance Setting Guide

This is **NOT** the pressure setting guide for **CODMAN CERTAS Plus Programmable Valves** (released in 2015).  
This guide applies only to CODMAN CERTAS Programmable Valves implanted before May 2013.



**Settings 1, 2, 3** ±20 mmH<sub>2</sub>O

**Setting 4** ±25 mmH<sub>2</sub>O

**Settings 5,6** ±35 mmH<sub>2</sub>O

**Settings 7** ±50 mmH<sub>2</sub>O

**Setting 8** is intended to limit flow. Average pressure > 400 mmH<sub>2</sub>O for flow rates 5 to 50ml/hr

## CODMAN NEURO



Johnson & Johnson Medical Ltd. Pinewood Campus, Nine Mile Ride,  
Wokingham, RG40 3EW United Kingdom  
[www.depuyshnthes.com](http://www.depuyshnthes.com)

Codman Neuro is a trading division of Johnson & Johnson Medical Ltd.  
**Registered Office:** PO Box 1988, Simpson Parkway, Livingston,  
West Lothian, EH54 0AB, Registered in Scotland, No.SC 132162